



## SPECIALIST

---

As we manufacture a soft acrylic for the production of hearing aids, we made this section for all of our non-dental customers.

### IN THIS SECTION:

Acrylic	454
Duplicating	457
Equipment	459
Miscellaneous	463
Separator	463

## Acrylic



### CENTRI™ OPACRYL

Proudly made by WHW, this excellent quality cold cure orthodontic resin does not run or slump and offers exceptional value for money. The material was developed primarily for the salt and pepper technique and has great strength and clarity.

Code	Description	Size
01OP1	Opacryl Powder	1kg
01OP3	Opacryl Powder	3kg
01OL1R	Opacryl Liquid Rose	1ltr
01OL1C	Opacryl Liquid Clear	1ltr
01OL50R	Opacryl Liquid Rose	500ml
01OL50C	Opacryl Liquid Clear	500ml
01OSPR	Opacryl Powder & Rose Liquid Kit	1kg & 500ml
01OSPC	Opacryl Powder & Clear Liquid Kit	1kg & 500ml

### CENTRI™ OPACRYL COLOUR LIQUIDS

Extensive range of colours that can be used as is or with our Centri™ Opacryl clear liquid to create a more translucent effect.

Code	Description	Size
01OCLBK	Opacryl Colour Black	125ml
01OCLBL	Opacryl Colour Blue	125ml
01OCLGR	Opacryl Colour Green	125ml
01OCLOR	Opacryl Colour Orange	125ml
01OCLRD	Opacryl Colour Red	125ml
01OCLRB	Opacryl Colour Rubine	125ml
01OCLTQ	Opacryl Colour Turquoise	125ml
01OCLVT	Opacryl Colour Violet	125ml
01OCLYW	Opacryl Colour Yellow	125ml

## Acrylic

### CENTRI™ SOFT ACRYLIC

Proudly made by WHW, this excellent heat cure soft acrylic has been manufactured for the production of hearing aids and other non-dental appliances.

Code	Description	Size
01SAAMP	Powder	3kg
01SAAM50	Powder	45kg
01SAAM1	Liquid	1ltr



### ORTHOCRYL® POWDER

Clear cold curing orthodontic resin especially formulated for the salt and pepper technique.

Code: 38160-112-00 Size: 1kg



### ORTHOCRYL® EQ POWDER

Clear cold curing orthodontic resin especially formulated for the doughing technique.

Code: 38160-300-00 Size: 1kg

### ORTHOCRYL® LIQUID

Standard cold curing liquids for both the Orthocryl® and Orthocryl® EQ powders.

Code	Description	Size
38161-100-00	Orthocryl® Liquid Clear	500ml
38161-300-00	Orthocryl® Liquid Transparent Pink	500ml
38161-150-00	Orthocryl® Liquid Clear	1ltr
38161-350-00	Orthocryl® Liquid Transparent Pink	1ltr



## Acrylic

### ORTHOCRYL® COLOUR LIQUIDS

Standard and neon cold curing liquids for both the Orthocryl® and Orthocryl® EQ powders.



Code	Description	Size
38161-127-00	Orthocryl® Liquid Red	250ml
38161-128-00	Orthocryl® Liquid Green	250ml
38161-129-00	Orthocryl® Liquid Blue	250ml
38161-130-00	Orthocryl® Liquid Yellow	250ml
38161-131-00	Orthocryl® Liquid Turquoise	250ml
38161-132-00	Orthocryl® Liquid Pink	250ml
38161-133-00	Orthocryl® Liquid Emerald Green	250ml
38161-134-00	Orthocryl® Liquid Violet	250ml
38161-135-00	Orthocryl® Liquid Neon Yellow	250ml
38161-136-00	Orthocryl® Liquid Neon Orange	250ml
38161-137-00	Orthocryl® Liquid Neon Pink	250ml
38161-138-00	Orthocryl® Liquid Neon Green	250ml



### ORTHOCRYL® COLOUR CONCENTRATE

Colour effects for tinting clear Orthocryl® liquid.

Code	Description	Size
38161-620-00	Orthocryl® Concentrate Red	100ml
38161-622-00	Orthocryl® Concentrate Blue	100ml



### ORTHOCRYL® BLACK & WHITE

The Orthocryl® black & white powder technique opens up a world of new ideas and creative possibilities for orthodontic appliances.

Code	Description	Size
38160-011-00	Orthocryl® Powder Black	200g
38160-012-00	Orthocryl® Powder White	200g
38160-013-00	Orthocryl® Powder Black	1kg
38160-014-00	Orthocryl® Powder White	1kg
38161-100-00	Orthocryl® Liquid Clear	500ml
38161-150-00	Orthocryl® Liquid Clear	1ltr



## Acrylic

### ORTHOCRYL® DISCO GLITTER

Incredibly good value, this glitter is added to the powder and processed as normal. Add the desired amount to achieve the required lustre.

Code	Description	Size
38161-605-00	Iris	50g
38161-606-00	Red	50g
38161-607-00	Blue	50g
38161-608-00	Silver	50g
38161-609-00	Gold	50g
38161-611-00	Mother of Pearl	50g



## Duplicating

### PLASTIC DUPLICATING FLASK

Suitable for use with Agar-Agar/duplicating gels.

Code: 51PDF

Size: 65 x 100 x 90mm



### SILICONE DUPLICATING FLASK

Adjustable duplicating flask for use with silicone.

Code: 01ZDF

Size: 45 x 90 x 75mm



### DENTSPLY SIRONA POLYFLEX®

Highest quality universal duplicating gel with properties comparable to silicone. Tear resistant with incredibly high definition.

Code: 95H5525300195

Size: 7.5kg



## Duplicating

### DUBLITOP®

Exceptional quality duplicating gel with extremely high tear resistance and detail reproduction. Typical amount of reuses - 25. Light blue - not suitable for plaster models.

Code: 38108-610-00

Size: 5kg



### DUPLIKAT

Top quality duplicating gel with very high tear resistance and detail reproduction. Typical amount of reuses - 20. Yellow - not suitable for plaster models.

Code: 38108-600-00

Size: 6kg



### DUBLIFORM

Time-tested duplicating gel with high tear resistance and detail reproduction. Typical amount of reuses - 15. Orange - suitable for use in the microwave but not suitable for plaster models.

Code: 38108-500-00

Size: 6kg



### DUBLIPLAST®

Universal duplicating gel with high tear resistance and detail reproduction. Typical amount of reuses - 15. Violet - suitable for plaster models

Code: 38165-500-00

Size: 5kg



## Duplicating



### CENTRI™ DUPLIKA 12

Ever popular purple duplicating silicone with high tear resistance and incredible definition.

Code	Size
68CD12S	1kg + 1 kg
68CD12L	5kg + 5kg

### CENTRI™ DUPLIKA HARD

Pink duplicating silicone with a denser matrix providing even higher tear resistance and definition.

Code	Size
68CDHS	1kg + 1kg
68CDHL	5kg + 5kg

## Equipment

### OMEC POLISHING LATHE

Two speed (1400-2800 rpm), mountable, heavy duty lathe with 350W motor.


Code	Size	Dimensions
60OPL	10kg	210 x 360 x 210mm



## Equipment

### IEU BENCH TOP POLISHING UNIT

Constructed in galvanised steel with a white, tough plastic coating for added durability. The powerful dust extraction system is driven by two centrifugal fans (310 cfm each), which power the troughs independently. These troughs are manufactured in polypropylene and are fitted with a clear adjustable plastic visor. All particles are extracted through our efficient pre and main filtration system, these are easy and cost effective to replace. The top of the unit is hinged and offers easy access to the filters and fluorescent lighting system.



Code	Size	Dimensions
01IEU	40kg	600 x 750 x 515mm

### POLISHING UNIT CONSUMABLES

Code	Description	Pack
01IEUFB	Dust Bags	2pcs
01IEUFBE	Exhaust Filters	2pcs
01IEUFP	Filter Pads	2pcs

### DREVE POLYMAX RANGE

Suitable for hot and cold polymerisation.

- ▶ 90 minute timer
- ▶ Adjustable air pressure >5 bar (Polymax 5 >6 bar)
- ▶ Automatic safety sealing system
- ▶ Maintenance free
- ▶ Drainage valve

### DREVE POLYMAX 1 - 95

Maximum Temperature	95°C
Pressure Pot	160ø x 110mm
Capacity	1 x brass flask or approximately 2 x pourable flasks

Code	Description	Size	Dimensions
D3429	Red	10.2kg	315 x 295 x 225mm
D34291	Silver	10.2kg	315 x 295 x 225mm



## Equipment

### DREVE POLYMAX 1 - 120

Maximum Temperature	120°C
Pressure Pot	160ø x 110mm
Capacity	1 x brass flask or approximately 2 x pourable flasks

Code	Description	Size	Dimensions
D3428	Red	10.2kg	315 x 295 x 225mm
D34281	Silver	10.2kg	315 x 295 x 225mm



### DREVE POLYMAX 3

Maximum Temperature	95°C
Pressure Pot	160ø x 300mm
Capacity	3 x brass flasks or approximately 6 x pourable flasks

Code	Description	Size	Dimensions
D3422	Red	15.1kg	315 x 295 x 440mm
D34221	Silver	15.1kg	315 x 295 x 440mm



### DREVE POLYMAX 5

Maximum Temperature	95°C
Pressure Pot	220ø x 174mm
Capacity	1 x large articulator, 3 x brass flasks or approximately 6 x pourable flasks

Code	Description	Size	Dimensions
D3423	Red	16.9kg	385 x 345 x 310mm
D34231	Silver	16.9kg	385 x 345 x 310mm



### BOSWORTH® PRESSURE POT

Hydraulic water press (1 litre capacity) used for creating porosity free relines, repairs, splints or any self-curing acrylic with just 10-15 Psi.

Code	Size	Dimensions
71C092135	2.8kg	127ø x 88mm (bowl)





## Equipment

### POLYCLAV® PRESSURE VESSEL

With a capacity of 4 litres, the easy to use Polyclav® is suitable for the polymerisation of all cold curing devices.

Code	Size	Dimensions
38070-000-00	3.2kg	200ø x 180mm (bowl)



### PRESSURE PAIL

Internal pressure can be generated either hydraulically using water or pneumatically using air from a compressor. The Pressure Pail is large enough to hold a standard relining jig or smaller articulator. The easy crank handle, replaceable O-ring and all-aluminium construction make this pot excellent for curing all types of cold cured acrylics. Great for repairs, relines, splint production, retainers and more. **Airline adaptor sold separately.**

Code	Size	Dimensions
399622760	2.3kg	152ø x 152mm (bowl)



Code	Description	Pack
399622780	Airline adaptor	1pc

### TRAUS B110

Universal brushless micromotor with foot, knee and table top control. Featuring cruise control, reverse, digital display and rpm control dial.



Code	Size	Dimensions
TSB110	3kg	235 x 90 x 200mm

### MARATHON ECO 450

Economic foot controlled, brushed micromotor.



Code	Size	Dimensions
63ECO450	2.3kg	185 x 100 x 140mm

## Miscellaneous

### AMANN GIRRBACH ARBOR BANDS

High quality and long-lasting.

Code	Description	Pack
AG815310	Coarse (80 grit)	50pcs
AG815330	Medium (120 grit)	50pcs
AG815340	Fine (240 grit)	50pcs
AG815300	Mandrel	1pc



### ARBOR BANDS

Cost effective and efficient, in a range of our grits and sizes.

Code	Description	Pack
07ABC12	Coarse 1/2"	100pcs
07ABC14	Coarse 1/4"	100pcs
07ABF12	Fine 1/2"	100pcs
07ABF14	Fine 1/4"	100pcs
07ABM12	Medium 1/2"	100pcs
07ABM14	Medium 1/4"	100pcs
07MAB12	Mandrel 1/2"	1pc
07MAB14	Mandrel 1/4"	1pc



## Separator

### CENTRI™ PCS

Plaster coating solution.

Code	Size
01PLC6	650ml
01PLC4	5ltr
01PLC25	25ltr

